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1. McMURRAY JJ, et al., ESC Committee for Practice Guidelines. ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure 2012: The Task Force for the Diagnosis and Treatment of Acute and Chronic Heart Failure 2012 of the European Society of Cardiology. Developed in collaboration with the Heart Failure Association (HFA) of the ESC. Eur Heart J. 2012;33:1787-847.
2. HILDEBRANDT P, COLLINSON PO, et al. Age-dependent values of N-terminal pro-B-type natriuretic peptide are superior to a single cut-point for ruling out suspected systolic dysfunction in primary care. Eur Heart J. 2010;31:1881-9.
3. DANIELS LB, MAISEL AS. Natriuretic peptides. J Am Coll Cardiol. 2007;50:2357-68.
4. VAN KIMMENADE RR, JANUZZI JL JR, et al. Renal clearance of B-type natriuretic peptide and amino terminal pro-B-type natriuretic peptide a mechanistic study in hypertensive subjects. J Am Coll Cardiol. 2009;53:884-90.
5. KIM HN, JANUZZI JL Jr. Natriuretic peptide testing in heart failure. Circulation. 2011;123:2015-9.
6. EMDIN M, PASSINO C, et al. Comparison of brain natriuretic peptide (BNP) and amino-terminal ProBNP for early diagnosis of heart failure. Clin Chem. 2007;53:1289-97.
7. O'DONOGHUE M, CHEN A, et al. The effects of ejection fraction on N-terminal ProBNP and BNP levels in patients with acute CHF: analysis from the ProBNP Investigation of Dyspnea in the Emergency Department (PRIDE) study. J Card Fail. 2005;11(Suppl.):S9-14.
8. HUNT SA, et al. 2009 focused update incorporated into the ACC/AHA 2005

Guidelines for the Diagnosis and Management of Heart Failure in Adults: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines: developed in collaboration with the International Society for Heart and Lung Transplantation. Circulation. 2009;119:e391-479.

9. MOE GW, HOWLETT J, et al.; Canadian Multicenter Improved Management of Patients With Congestive Heart Failure (IMPROVE-CHF) Study Investigators. N-terminal pro-B-type natriuretic peptide testing improves the management of patients with suspected acute heart failure: primary results of the Canadian prospective randomized multicenter IMPROVE-CHF study. Circulation. 2007;115:3103-10.
10. RUTTEN JH, STEYERBERG EW, et al.. N-terminal pro-brain natriuretic peptide testing in the emergency department: beneficial effects on hospitalization, costs, and outcome. Am Heart J. 2008;156:71-7.
11. JANUZZI JL, VAN KIMMENADE R, et al. NT-proBNP testing for diagnosis and short-term prognosis in acute destabilized heart failure: an international pooled analysis of 1256 patients: the International Collaborative of NT-proBNP Study. Eur Heart J. 2006;27:330-7.
12. BAGGISH AL, VAN KIMMENADE RR, JANUZZI JL JR. The differential diagnosis of an elevated amino-terminal pro-B-type natriuretic peptide level. Am J Cardiol. 2008;101(Suppl.):43A-8A. 1. Guidelines for the diagnosis and treatment of Chronic Heart Failure (update 2005). The European Society of Cardiology
13. BAGGISH AL, VAN KIMMENADE RR, JANUZZI JL Jr. Amino-terminal pro-B-type natriuretic peptide testing and prognosis in patients with acute dyspnea, including those with acute heart failure. Am J Cardiol. 2008;101(Suppl.):49A- 55A.
14. MASSON S, LATINI R. Amino-terminal pro-B-type natriuretic peptides and prognosis in chronic heart failure. Am J Cardiol. 2008;101 (Suppl.):56A-60A.
15. OMLAND T, DE LEMOS JA. Amino-terminal pro-B-type natriuretic peptides in stable and unstable ischemic heart disease. Am J Cardiol. 2008;101(Suppl.):61A-6.
16. WEBER M, BAZZINO O, et al. N-terminal B-type natriuretic peptide assessment provides incremental prognostic information in patients with acute coronary syndromes and normal troponin T values upon admission. J Am Coll Cardiol. 2008;51:1188-95.
17. BETTENCOURT P, JANUZZI JL Jr. Amino-terminal pro-B-type natriuretic peptide testing for inpatient monitoring and treatment guidance of acute destabilized heart failure. Am J Cardiol. 2008;101(Suppl.):67A-71A.
18. PORAPAKKHAM P, PORAPAKKHAM P, et al. B-type natriuretic peptide-guided heart failure therapy: A metaanalysis. Arch Intern Med. 2010;170:507-14.
19. PASSING H, BABLOK W, et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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