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
---

**Run Information**

Sample ID	Example Report	Run Date	12 Jul 2016 12:00 AM
Protocol	Sputum v3.3	Serial No.	01234567
Pouch Type	Pneumo v2.0	Lot No.	012345
Controls	Passed	Operator	Anonymous
Run Status	Completed	Instrument	FA0000

---

**Detection Summary**


**Bacteria**

	Bin (copies/mL)	Bin (copies/mL)			
		$10^4$	$10^5$	$10^6$	$\geq 10^7$
Detected: <b>3a</b>	$\geq 10^7$				
	$10^6$				
	$10^4$				

 Note: Detection of bacterial nucleic acid may be indicative of colonizing or normal respiratory flora and may not indicate the causative agent of pneumonia. Semi-quantitative Bin (copies/mL) results generated by the FilmArray Pneumonia Panel are not equivalent to CFU/mL and do not consistently correlate with the quantity of bacterial analytes compared to CFU/mL. For specimens with multiple bacteria detected, the relative abundance of nucleic acids (copies/mL) may not correlate with the relative abundance of bacteria as determined by culture (CFU/mL). Clinical correlation is advised to determine significance of semi-quantitative Bin (copies/mL) for clinical management.

---

**Antimicrobial Resistance Genes**

Detected:  CTX-M

 Note: Antimicrobial resistance can occur via multiple mechanisms. A Not Detected result for a genetic marker of antimicrobial resistance does not indicate susceptibility to associated antimicrobial drugs or drug classes. A Detected result for a genetic marker of antimicrobial resistance cannot be definitively linked to the microorganism(s) detected. Culture is required to obtain isolates for antimicrobial susceptibility testing and FilmArray Pneumonia Panel results should be used in conjunction with culture results for the determination of susceptibility or resistance.

---

**Atypical Bacteria**

Detected: None **4**

---



**Viruses**

Detected:  Influenza A

---

**Result Summary**


**Bacteria**

	Bin (copies/mL)	Bin (copies/mL)			
		$10^4$	$10^5$	$10^6$	$\geq 10^7$
Not Detected	<i>Acinetobacter calcoaceticus-baumannii</i> complex				
Not Detected	<i>Enterobacter cloacae</i> complex				
Not Detected	<i>Escherichia coli</i>				
	$10^4$				
Not Detected	<i>Klebsiella aerogenes</i>				
Not Detected	<i>Klebsiella oxytoca</i>				
	$\geq 10^7$				
Not Detected	<i>Moraxella catarrhalis</i>				
Not Detected	<i>Proteus</i> spp.				
Not Detected	<i>Pseudomonas aeruginosa</i>				
Not Detected	<i>Serratia marcescens</i>				
Not Detected	<i>Staphylococcus aureus</i>				
Not Detected	<i>Streptococcus agalactiae</i>				
Not Detected	<i>Streptococcus pneumoniae</i>				
	$10^6$				

 Note: Detection of bacterial nucleic acid may be indicative of colonizing or normal respiratory flora and may not indicate the causative agent of pneumonia. Semi-quantitative Bin (copies/mL) results generated by the FilmArray Pneumonia Panel are not equivalent to CFU/mL and do not consistently correlate with the quantity of bacterial analytes compared to CFU/mL. For specimens with multiple bacteria detected, the relative abundance of nucleic acids (copies/mL) may not correlate with the relative abundance of bacteria as determined by culture (CFU/mL). Clinical correlation is advised to determine significance of semi-quantitative Bin (copies/mL) for clinical management.


---

**Antimicrobial Resistance Genes**

 Detected CTX-M

Not Detected IMP

Not Detected KPC

 N/A *mecA/C* and MREJ

Not Detected NDM

Not Detected OXA-48-like

Not Detected VIM

 Note: Antimicrobial resistance can occur via multiple mechanisms. A Not Detected result for a genetic marker of antimicrobial resistance does not indicate susceptibility to associated antimicrobial drugs or drug classes. A Detected result for a genetic marker of antimicrobial resistance cannot be definitively linked to the microorganism(s) detected. Culture is required to obtain isolates for antimicrobial susceptibility testing and FilmArray Pneumonia Panel results should be used in conjunction with culture results for the determination of susceptibility or resistance.

---

**Atypical Bacteria**

Not Detected *Chlamydia pneumoniae*

Not Detected *Legionella pneumophila*

Not Detected *Mycoplasma pneumoniae*

---

**Viruses**

Not Detected Adenovirus

Not Detected Coronavirus

Not Detected Human Metapneumovirus

Not Detected Human Rhinovirus/Enterovirus

 Detected Influenza A

Not Detected Influenza B

Not Detected Parainfluenza Virus

Not Detected Respiratory Syncytial Virus