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1. Balows A, Hausler WJ, Herrmann KL, Isenberg HD, Shadomy HJ, editors. *Manual of Clinical Microbiology*, 5th ed. American Society for Microbiology, Washington, D.C. 1991.
2. Balows A, Truper HG, Dworkin M, Harder W, and Schleifer K-H, editors. *The Prokaryotes - a Handbook on the Biology of Bacteria: Exophysiology, Isolation, Identification, Applications*, 2nd ed., Volume II. Springer-Verlag, New York 1992.

3. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
4. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C 263a. PL 100-578. 1988.
5. Difco Manual Dehydrated Culture Media and Reagents for Microbiology. 10th ed. 1984.
6. Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. *Bergey's Manual of Determinative Bacteriology*, 9th ed. Williams and Wilkins, Baltimore 1994.
7. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. *Color Atlas and Textbook of Diagnostic Microbiology* 4th ed. Lippincott, Philadelphia, PA 1992.
8. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. *Color Atlas and Textbook of Diagnostic Microbiology* 5th ed. Lippincott, Philadelphia, PA. 1997.
9. Krieg NR, Holt JG, editors. *Bergey's Manual of Systematic Bacteriology*, 9th ed. Williams and Wilkins, Baltimore 1984.
10. Manafi M, Kneifel W, Bascomb S. Fluorogenic and chromogenic substrates used in bacterial diagnostics. *Microbiol.* 1991; Rev. 55: 335-348.
11. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Yolkner RH, editors. *Manual of Clinical Microbiology*, 7th ed. American Society for Microbiology, Washington, D.C. 1999.
12. Murray PR, Baron EJ, Jorgensen JH, Pfaller MA and Yolkner RH, editors. *Manual of Clinical Microbiology*, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
13. National Committee for Clinical Laboratory Standards, M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue* - Approved Guideline, 1997.
14. Nørskov-Lauritsen N, Kilian M. 2006. Reclassification of *Actinobacillus actinomycetemcomitans*, *Haemophilus aphrophilus*, *Haemophilus paraphrophilus*, and *Haemophilus segnis* as *Aggregatibacter actinomycetemcomitans* gen. nov., comb. nov., *Aggregatibacter aphrophilus* comb. nov., and *Aggregatibacter segnis* comb. nov., and emended description of *Aggregatibacter aphrophilus* to include V factor-dependent and V factor-independent isolates. *IJSEM.* 2006. 56:2135-2146.
15. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, *Biosafety in Microbiological and Biomedical Laboratories*, 1988.
16. Weyant RS, Moss CW, Weaver RE, Hollis DG, Jordon JG, Cook EC, and Daneshvar MI. Identification of Unusual Pathogenic and Gram-Negative Aerobic and Facultatively Anaerobic Bacteria 2nd ed. Williams & Wilkins, Philadelphia, PA 1996.

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Product covered by one or more of U.S. Patent Numbers D 414,272 ; D 437,797 ; 5,609,828 ; 5,746,980 ; 5,804,437 ; 5,869,005 ; 5,932,177 ; 5,951,952 ; 6,267,929 ; 6,309,890 and 6,340,573 ; and Foreign Counterparts.
Other Patents Pending.

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