フラザイオン(FRASER-T)*For microbiological control only***FRASER broth (FRASER-T)**

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1. NF EN ISO 11290-1/A1 (February 2005) - Microbiology of food and animal feeding stuffs - Horizontal method for the detection and enumeration of *Listeria monocytogenes* - Part 1 : detection method - Amendment 1 : modification of the isolation media, of the haemolysis test and inclusion of precision data.
2. FRASER J.A., SPERBER W.H. - Rapid Detection of *Listeria* spp. in Food and Environmental Samples by Esculin Hydrolysis - *J. of Food protection*, 1988, vol. 51, p.762-765.
3. RODRIGUEZ L, FERNANDEZ G., GARAYZABAL J. et al. - New Methodology for the Isolation of *Listeria* Microorganisms from Heavily Contaminated Environments. - *Applied and Environmental Microbiology*, 1984, vol. 47, p. 1188-1190.
4. HOLBROOK R., ANDERSON J.M., BRIGGS T.A. et al. - Faster detection of *Listeria* in food using rapid immunoassay following culture - 3rd World Congress foodborne infections and intoxications, (16-19 June 1992), Berlin, p.1208-1210.
5. XP CEN ISO/TS 11133-2 (2004) - Microbiology of food and animal feeding stuffs. Guidelines on preparation and production of culture media. Part 2: Practical guidelines on performance testing of culture media.

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