chromID™ CARBA (CARB)

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

1. CDC - Guidance for Control of infections with carbapenem-resistant or carbapenemase-producing *Enterobacteriaceae* in acute care facilities - MMWR, 2009, 58(10): 256-260.
2. PERRY J. D., NAQVI Sakeenah Hussain, MIRZA Irfan Ali et al. - Prevalence of faecal carriage of *Enterobacteriaceae* with NDM-1 carbapenemase at military hospitals in Pakistan, and evaluation of two chromogenic media. - Journal of Antimicrobial Chemotherapy, 2011, ISSN:1460-2091.
3. WALSH T.R. - Emerging carbapenemases: a global perspective. International Journal of Antimicrobial Agents, 2010, 36S3: 8-14.
4. NORDMANN P., CUZON G., NAAS T. - The real threat of *Klebsiella pneumoniae* carbapenemase-producing bacteria. Lancet Infectious Disease, 2009, 9: 228-236.
5. ORENGA S., JAMES A.L., PERRY J.D., PINCUS D.H. (2009). Enzymatic substrates in microbiology. Journal of Microbiological Methods, 79: 139-155.
6. Laboratory protocol for detection of carbapenem-resistant or carbapenemase-producing *Klebsiella* spp. and *E. coli* from rectal swabs
-www.cdc.gov/HAI/pdfs/labSettings/Klebsiella_or_Ecoli.pdf
7. VRIONI G., DANIIL I., VOULGARI E., RANELLOU K., KOUMAKI V., GHIRARDI S., KIMOULI M., ZAMBARDI G. and TSAKRIS A. - Comparative evaluation of a prototype medium (chromID CARBA) for detecting carbapenemase-producing *Enterobacteriaceae* in surveillance rectal swabs. *J. Clin. Microbiol.*, June 2012, in press.
8. BEREKSI N., GIRAUD D., JOYEUX F. and al. - Evaluation of a chromogenic medium, chromID CARBA, for the detection of carbapenemase-producing *Enterobacteriaceae*. - Poster 1718 - London 2012 - 22nd ECCMID.
9. VRIONI G., DANIIL J., VOULGARI E. and al. - Evaluation of a novel chromogenic medium for detecting carbapenemase-producing *Enterobacteriaceae* in surveillance rectal swabs - Poster 1717 - London 2012 - 22nd ECCMID.
10. WILKINSON K., ARMES A., RAZA M. and PERRY J.D. - Evaluation of various culture media and procedures recommended for isolation of carbapenemase-producing *Enterobacteriaceae*. - Poster 2316 - London 2012 - 22nd ECCMID.

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Fax 33 (0)4 78 87 20 90

www.biomerieux.com

The logo is a registered and protected trademark of bioMérieux sa or one of its subsidiaries.