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Singer M., et al. JAMA. 2016;315(8):801-810.

<http://jama.jamanetwork.com/article.aspx?articleid=2492881>

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日本版敗血症診療ガイドライン2016 (The Japanese Clinical Practice Guidelines for Management of Sepsis and Septic Shock 2016 (J-SSCG2016))

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1. Singer M., *et al.* The Third International Consensus Definitions for Sepsis and Septic Shock (Sepsis-3). *JAMA*. 2016;315(8):801-810
2. Koneman E.W., *et al.*, Color Atlas and Textbook of Diagnostic Microbiology. Third Edition
3. Principles and procedures for Blood Cultures; Approved Guideline, CLSI document M47-A. Clinical and Laboratory Standards Institute (CLSI); Wayne, P.A. 2007
4. Adhikari N.K.J., Fowler R.A., Bhagwanjee S., Rubenfeld G.D., Critical care and the global burden of critical illness in adults. *Lancet* 2010;376:1339–1346
5. WSD fact sheet 2013/www.world-sepsis-day.org
6. Kumar, A, *et al.* Initiation of inappropriate antimicrobial therapy results in a fivefold reduction of survival in human septic shock. *Chest*. 2009 Nov;136(5):1237-48
7. European Society for Clinical Microbiology and Infectious Diseases, European Manual for Clinical Microbiology, 1st Edition, 2012
8. Garey KW., Rege M., Manjunath P. Pai, Mingo DE., Suda KJ., Turpin RS., Bearden DT. Time to Initiation of Fluconazole Therapy Impacts Mortality in Patients with Candidemia: A Multi-Institutional Study. *Clin Infect Dis*. 2006;43(1):25-31
9. Khatib R., Saeed S., Sharma M., Riederer K., Fakh MG., Johnson LB. Impact of initial antibiotic choice and delayed appropriate treatment on the outcome of *Staphylococcus aureus* bacteremia. *Eur J Clin Microbiol Infect Dis*. 2006;25(3):181-5
10. Kollef MH, Sherman G, Ward S, Fraser VJ. Inadequate antimicrobial treatment of infections: a risk factor for hospital mortality among critically ill patients. *Chest*. 1999;115(2):462-74
11. Harbarth S., Garbino J., Pugin J., Romand J.A., Lew D., Pittet D. Inappropriate initial antimicrobial therapy and its effect on survival in a clinical trial of immunomodulating therapy for severe sepsis. *Am J Med*. 2003;115(7):529–535
12. Lodise T.P., McKinnon P.S., Swiderski L., Rybak M.J. Outcomes Analysis of Delayed Antibiotic Treatment for Hospital-Acquired *Staphylococcus aureus* Bacteremia. *CID* 2003;36:1419-1423
13. Kang C.I., Kim S.H., Kim H.B., Park S.W., Choe Y.J., Oh M.D., Kim E.C., Choe K.W. *Pseudomonas aeruginosa* bacteremia: risk factors for mortality and influence of delayed receipt of effective antimicrobial therapy on clinical outcome. *Clin Infect Dis*. 2003;37(6):745-51
14. Forrest G.N., Mankes K., Jabra-Rizk M.A., Weekes E., Johnson J.K., Lincalis D.P., Venezia R.A. Peptide Nucleic Acid Fluorescence In Situ Hybridization-Based Identification of *Candida albicans* and Its Impact on Mortality and Antifungal Therapy Costs. *J Clin Microbiol*. 2006 Sep; 44(9): 3381–3383
15. Kumar A, *et al.* Duration of hypotension before initiation of effective antimicrobial therapy is the critical determinant of survival in human septic shock. *Crit Care Med*. 2006;34(6):1589-96
16. Baron, E.J., M.P. Weinstein, W.M. Dunne, Jr., P. Yagupsky, D.F. Welch, and D.M. Wilson. Cumitech 1C, Blood Cultures IV. Coordinating ed., E.J. Baron. ASM Press, Washington, D.C. 2005
17. Mermel L.A., Maki D.G. Detection of bacteremia in adults : consequences of culturing an inadequate volume of blood. *Ann Intern Med*. 1993;119:270-272
18. Bouza E, Sousa D, Rodriguez-Cr eixems M, Lechuz JG, Mu oz P. Is the volume of blood cultured still a significant factor in the diagnosis of bloodstream infections? *J Clin Microbiol*. 2007 45:2765-9
19. Cockerill FR III, Wilson J.W., Vetter E.A., *et al.* Optimal testing parameters for blood cultures. *Clin Infect Dis*. 2004 ;38 :1724-1730
20. Kellogg J.A., Manzella J.P., Bankert D.A. Frequency of low-level bacteremia in children from birth to fifteen years of age. *J Clin Microbiol*. 2000;38:2181-2185
21. Freedman S.B., Roosevelt G.E. Utility of anaerobic blood cultures in a pediatric emergency department. *Pediatr Emerg Care*. 2004;20(7):433-6
22. Lee A., Weinstein MP., Mirrett S., Reller LB. Detection of Bloodstream Infections in Adults: How Many Blood Cultures Are Needed? *J Clin Microbiol* 2007;45:3546-3548
23. Lee DH., Kim S.C., Bae IG., Koh EH., Kim S., Clinical Evaluation of BacT/ALERT FA Plus and FN Plus Bottles Compared with Standard Bottles , *J. Clin. Microbiol*. 2013; 51(12): 4150-4155
24. Amarsy-Guerle R., Mougari F., Jacquier H., Oliari J., Benmansour H., Riahi J., Ber ot B., Raskine L., Cambau E., High medical impact of implementing the new polymeric bead-based BacT/ALERT[ ] FA Plus and FN Plus blood culture bottles in standard care, *Eur J. Clin. Microbiol Dis*. 2015;34(5):1031-1037
25. Kirn T.J., Mirrett S., Reller L.B., Weinstein M.P., Controlled Clinical Comparison of BacT/ALERT FA Plus and FN Plus Blood Culture Media with BacT/ALERT FA and FN Blood Culture Media, *J. Clin. Microbiol*. 2014; 52(3): 839-843
26. Doern C., Mirrett S., Halstead D., Abid J., Okada P., Reller L.B. Controlled Clinical Comparison of New Pediatric Medium with Adsorbent Polymeric Beads (PF Plus) versus Charcoal-Containing PF Medium in the BacT/ALERT Blood Culture System, *J. Clin. Microbiol*. 2014; 52(6): 1898-1900
27. Riley J.A., Heiter B.J., Bourbeau P.P. Comparison of recovery of blood culture isolates from two BacT/ALERT FAN aerobic blood culture bottles with recovery from one FAN aerobic bottle and one FAN anaerobic bottle. *J.Clin. Microbiol*. 2003;41:213-217
28. Weinstein, M. P., M. L. Towns, S. M. Quartey, S. Mirrett, L. G. Reimer, G. Parmigiani, and L. B. Reller. The clinical significance of positive blood cultures in the 1990s; a prospective comprehensive evaluation of the microbiology, epidemiology, and outcome of bacteremia and fungemia in adults. *Clin. Infect. Dis*. 1997;24:584–602
29. UK Department of Health: Taking Blood Cultures – A summary of best practice. 2007
30. Weinstein M.P. Current blood culture methods and systems: clinical concepts, technology, and interpretation of results. *Clin Infect Dis*. 1996 ;23 :40-46
31. Everts R.J., Vinson E.N., Adholla P.O., Reller L.B. Contamination of catheter-drawn blood cultures. *J. Clin Microbiol*. 2001;39:3393-3394
32. Cornaglia G., *et al.* European Manual of Microbiology. ESCMID-SFM 2012

33. Kirm T.J., Weinstein M.P. Update on blood cultures: how to obtain, process, report, and interpret. *Clin Microbiol Infect.* 2013;19(6):513–520
34. Wilson M.L., Mirrett S., Reller L.B. *et al.* Recovery of clinically important microorganisms from the BacT/ALERT blood culture system does not require testing for 7 days. *Diagn Microbiol Infect Dis.* 1993;16:31-34
35. Bourbeau PP., Foltzer M. Routine incubation of BacT/ALERT FA and FN blood culture bottles for more than 3 days may not be necessary. *J Clin Microbiol.* 2005;43:2506-2509
36. *Clinical Infectious Disease.* Edited by David Schlossberg. Cambridge University Press, 2015
37. Keri K. Hall and Jason A. Lyman, Updated Review of Blood Culture Contamination, *Clin. Microbiol. Rev.* 2006, 19(4):788
38. Dunne W.M. Jr., Nolte F.S., Wilson M.L. *Cumitech 1B, Blood Cultures III.* Coordinating ed., Hindler J.A. ASM Press. Washington, D.C. 1997
39. Hall, K.K. and J.A. Lyman. Updated review of blood culture contamination. *Clinical Microbiology Reviews.* 2006;19:788-802
40. Bamber, A.I., J. G. Cuniffe, D. Nayar, R. Ganguly and E. Falconer. The effectiveness of introducing blood culture collection packs to reduce contamination. *British Journal of Biomedical Science.* 2009;66(1):1-9.
41. Gander, R. M., L. Byrd, M. DeCrescenzo, S. Hirany and M. Bowen, J. Baughman. Impact of blood cultures drawn by phlebotomy on contamination rates and health care costs in a hospital emergency department. *J. Clin. Microbiol.* 2009;47:1021-1024
42. Richter S.S., Beekman S.E., Croco D.J., Koontz R.P., Pfaller M.A., Doern G.V. Minimizing the workup of blood culture contaminants: implementation and evaluation of a laboratory-based algorithm. *J Clin Microbiol.* 2002;40:2437-2444
43. Towns M.L., Reller L.B. Diagnostic methods: current best practices and guidelines for isolation of bacteria and fungi in infective endocarditis. *Infect Dis Clin N Am.* 2002;16:363-376
44. Osborn TM., Nguyen HB., Rivers EP. Emergency medicine and the surviving sepsis campaign: an international approach to managing severe sepsis and septic shock. *Ann Emerg Med* 2005;46:228-231
45. Rubenstein E., Lang R. Fungal endocarditis. *Eur Heart J.* 1995; 16(Suppl B):84-89
46. Ellis ME., Al-Abdely H., Standridge A., Greer W., Venture W. Fungal endocarditis: evidence in the world literature, 1965-1995. 2001;32:50-62
47. McLeod R., Remington JS. Fungal endocarditis. In: Rahimtoola SH *et al.*, eds. *Infective Endocarditis.* New York, NY: Gune & Stratton.1978:211-290
48. Ziegler R., Johnscher I., Martus P., Lendardt D., Just HM. Controlled Clinical Laboratory Comparison of Two Supplemented Aerobic and Anaerobic Media Used in Automated Blood Culture Systems To Detect Bloodstream Infections. *J Clin Microbiol.* 1998;36:657-661
49. Pohlman JK., Kirkley BA., Easley KA., Basille BA., Washington JA. Controlled Clinical Evaluation of BACTEC Plus Aerobic/F and BacT/ALERT Aerobic FAN Bottles for Detection of Bloodstream Infections. *J Clin Microbiol.* 1995;33:2856-2858
50. Baron EJ., Scott JD., Tompkins LS. Prolonged incubation and extensive subculturing do not increase recovery of clinically significant microorganisms from standard automated blood cultures. *Clin Infect Dis.* 2005;41:1677-1680
51. Beekmann SE., Diekema D.J., Chapin KC., Goern GV. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. *J Clin Microbiol.* 2003;41:3119-3125
52. Munson E., Diekema DJ., Beekmann SE., Chapin KC., Doern GV. Detection and treatment of bloodstream infection: laboratory reporting and antimicrobial management. *J Clin Microbiol.* 2003;41:495-497
53. Barenfanger J, Graham DR, Kolluri L, Sangwan G, Lawhorn J, Drake CA, Verhulst SJ, Peterson R, Moja LB, Ertmoed MM, Moja AB, Shevlin DW, Vautrain R, Callahan CD. Decreased Mortality Associated With Prompt Gram Staining of Blood Cultures, *Am J Clin Pathol* 2008;130:870-876
54. Timbrook T, Boger MS, Steed LL, Hurst JM, 2015. Unanticipated Multiplex PCR Identification of Polymicrobial Blood Culture Resulting in Earlier Isolation, Susceptibilities, and Optimization of Clinical Care. *J. Clin. Microbiol.* 2015;53(7):2371-3
55. Bauer KA, West JE, Balada Llasat J, Pancholi P, Stevenson KB, Goff DA. An Antimicrobial Stewardship Program’s Impact with Rapid Polymerase Chain Reaction Methicillin Resistant *Staphylococcus aureus* / *S. aureus* Blood Culture Test in Patients with *S. aureus* Bacteremia. *Clin. Infect. Dis.* 2010;51(9):1074-1080.
56. Dierkes C, Ehrenstein B, Siebig S, Linde H-J, Reischl U, Salzberger B. Clinical impact of a commercially available multiplex PCR system for rapid detection of pathogens in patients with presumed sepsis. *BMC Infect. Dis.* 2009;9(1):126
57. Weinstein MP. Blood Culture Contamination: Persisting Problems and Partial Progress. *J Clin Microbiol.* 2003;41:2275-2278
58. Weinstein MP., Towns ML, Quartey SM. *et al.* The clinical significance of positive blood cultures in the 1990s: a prospective comprehensive evaluation of the microbiology, epidemiology and outcome of bacteremia and fungemia in adults. *Clin Infect Dis.* 1997;24:584-602
59. *Applied Phlebotomy.* Dennis J. Ernst, Dennis J. Ernst (MT(ASCP)). Lippincott Williams & Wilkins, 2005
60. *Essentials Of Medical Laboratory Practice.* Constance L Lieseke, Elizabeth A Zeibig. F.A. Davis, 2012
61. Qamruddin A, Khanna N, Orr D. Peripheral blood culture contamination in adults and venipuncture technique: prospective cohort study. *J Clin Pathol.* 2008 61:509-13

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